



2.00 mg/L C from NIST KHP -- 40 mL

Item No. 12100

2.00 mg/L C NIST KHP, 40 mL NIST KHP CRMs.
Substitute for item 12101.

SPECIFICATIONS

Product Family	Process Water Consumables
Product Category	Total Organic Carbon
Item Type	CRM
Accreditations	ISO 17034
Analytical Technology	Total Organic Carbon Analyzer (TOC)
Lead Time (Business Days)	5
Storage	4 +/- 2C (Refrigerator)
Shelf Life from date of manufacture (months)	2
Components per Pkg	1
Container Type	40mL Vial w/septacap
Concentration	2.00 mg/L
Volume	40 mL
Weight	0.15 lbs
Instrument Manufacturer	Teledyne Tekmar, Sievers

×

Questions? We are happy to help.

Application	Pharmaceutical
-------------	----------------

Preservative	H3PO4
--------------	-------