



USP System Suitability Set - 40 mL set

Item No. 18000

USP/EP Bulk Water System Suitability Set Includes (1) Reagent Water Rw, (1) 0.5 mg/L C USP Sucrose, and (1) 0.5 mg/L USP 1,4-Benzoquinone in 40 mL vials. Replaces Sievers CSTD 31004-01. Substitute for item 18004.

SPECIFICATIONS

Product Family	Process Water Consumables
Product Category	System Suitability
Item Type	CRM
Accreditations	ISO 17034
Analytical Technology	Total Organic Carbon Analyzer (TOC)
Lead Time (Business Days)	1
Storage	4 +/- 2C (Refrigerator)
Shelf Life from date of manufacture (months)	2
Components per Pkg	3
Container Type	40mL Vial w/septacap
Concentration	0.500 mg/L
Weight	0.5 lbs
Instrument Manufacturer	Teledyne Tekmar, Sievers
Application	Pharmaceutical

Preservative

H3PO4

Instrument Model

Sievers 500, Sievers 800, Sievers 900
