



Labcare
de Colombia



PRODUCT CATALOGUE

VETERINARY RAPID TEST



Zhejiang Pushkang Biotechnology Co., Ltd



Veterinary Rapid Diagnostic Test Kit

Canine

No.	Product name	Abbr.	Format	Specimen	Specification
VE90001	Canine Distemper Virus Antigen Rapid Test	CDV Ag	Cassette	Secretions	10T
VE90002	Canine Distemper Virus Antibody Rapid Test	CDV Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90003	Canine Parvovirus Antigen Rapid Test	CPV Ag	Cassette	Feces	10T
VE90004	Canine Parvovirus Antibody Rapid Test	CPV Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90005	Canine Influenza Virus Antigen Rapid Test	CIV Ag	Cassette	Secretions	10T
VE90006	Canine Coronavirus Antigen Rapid Test	CCV Ag	Cassette	Feces	10T
VE90007	Canine Parainfluenza Virus Antigen Rapid Test	CPIV Ag	Cassette	Secretions	10T
VE90008	Canine Adenovirus II Antigen Rapid Test	CAV-II Ag	Cassette	Secretions	10T
VE90009	Canine Adenovirus I Antigen Rapid Test	CAV-I Ag	Cassette	Secretions	10T
VE90010	Canine C-reactive Protein Rapid Test	C. CRP	Cassette	Serum/Plasma	10T
VE90011	Toxoplasma Antibody Rapid Test	TOXO Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90012	Canine Heartworm Antigen Rapid Test	CHW Ag	Cassette	Whole Blood/Serum/Plasma	10T
VE90013	Canine Leishmania Antibody Rapid Test	C. LSH Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90014	Brucella Antibody Rapid Test	Bru Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90015	Canine Ehrlichia Antibody Rapid Test	C. EHR Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90016	Canine Leptospirosis Antibody Rapid Test	C. Lepto Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90017	Babesia Gibsoni Antibody Rapid Test	BG Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90018	Rabies Antigen Rapid Test	RAB Ag	Cassette	Secretions	10T
VE90019	Rabies Antibody Rapid Test	RAB Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90020	Canine Lyme Antibody Rapid Test	C. Lyme Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90021	Canine Relaxin Rapid Test	C. RLN	Cassette	Whole Blood/Serum/Plasma	10T
VE90022	Canine Giardia Antigen Rapid Test	C. GIA Ag	Cassette	Feces	10T
VE90023	Canine Anaplasma Antibody Rapid Test	C. ANA Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90024	Canine Parvovirus/Corona Virus Antigen Combo Rapid Test	CPV/CCV Ag	Cassette	Feces	10T
VE90025	Canine Distemper Virus/Parainfluenza Virus Antigen Combo Rapid Test	CDV/CPIV Ag	Cassette	Secretions	10T
VE90026	Canine Distemper Virus/Parvo Virus Antibody Combo Rapid Test	CDV/CPV Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90027	Canine Distemper Virus/Adeno Virus Antigen Combo Rapid Test	CDV/CAV Ag	Cassette	Secretions	10T
VE90028	Canine Parvovirus/Corona Virus/Giardia Antigen Combo Rapid Test	CPV/CCV/GIA Ag	Cassette	Feces	10T
VE90029	Canine Distemper Virus/Adeno Virus/Influenza Virus Antigen Combo Rapid Test	CDV/CAV/CIV Ag	Cassette	Secretions	10T
VE90030	Canine Ehrlichia/Anaplasma Antibody Combo Rapid Test	EHR/ANA Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90031	Canine Ehrlichia/Lyme/Anaplasma Antibody Combo Rapid Test	EHR/LYM/ANA Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90032	Canine Heartworm Antigen/Ehrlichia/Lyme/Anaplasma Antibody Combo Rapid Test	CHW Ag/EHR/LYM/ANA Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90033	Canine Ehrlichia/Babesia Gibsoni/Anaplasma Antibody Combo Rapid Test	EHR/BG/ANA Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90034	Canine Heartworm Antigen/Ehrlichia/Babesia Gibsoni/Anaplasma Antibody Combo Rapid Test	CHW Ag/EHR/BG/ANA Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90035	Canine Pancreatic Lipase Rapid Test	CPL	Cassette	Serum/Plasma	10T
VE90036	Canine Fecal Occult Blood Rapid Test	C. FOB	Cassette	Feces	10T

Feline

No.	Product name	Abbr.	Format	Specimen	Specification
VE90037	Feline Panleucopenia Virus Antigen Rapid Test	FPV Ag	Cassette	Feces	10T
VE90038	Feline Panleucopenia Virus Antibody Rapid Test	FPV Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90039	Feline Infectious Peritonitis Virus Antigen Rapid Test	FIPV Ag	Cassette	Feces/Ascites	10T
VE90040	Feline Infectious Peritonitis Virus Antibody Rapid Test	FIPV Ab	Cassette	Whole Blood/Serum/Plasma/Ascites	10T
VE90041	Feline Coronavirus Antigen Rapid Test	FCoV Ag	Cassette	Feces/Ascites	10T
VE90042	Feline Coronavirus Antibody Rapid Test	FCoV Ab	Cassette	Whole Blood/Serum/Plasma/Ascites	10T
VE90043	Feline Leukemia Virus Antigen Rapid Test	FeLV Ag	Cassette	Whole Blood/Serum/Plasma	10T
VE90044	Feline Immunodeficiency Virus Antibody Test	FIV Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90045	Feline Giardia Antigen Rapid Test	F. GIA Ag	Cassette	Feces	10T
VE90046	Feline Anaplasma Antibody Rapid Test	F. ANA Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90047	Toxoplasma Antibody Rapid Test	TOXO Ab	Cassette	Whole Blood/Serum/Plasma	10T
VE90048	Feline Herpesvirus Antigen Rapid Test	FHV Ag	Cassette	Secretions	10T
VE90049	Feline Calicivirus Antigen Rapid Test	FCV Ag	Cassette	Secretions	10T
VE90050	Feline Heartworm Antigen Rapid Test	FHW Ag	Cassette	Whole Blood/Serum/Plasma	10T
VE90051	Feline Immunodeficiency Virus Antibody/Leukemia Virus Antibody Combo Rapid Test	FIV Ab/FeLV Ag	Cassette	Whole Blood/Serum/Plasma	10T
VE90052	Feline Herpes Virus/Calici Virus Antigen Combo Rapid Test	FHV/FCV Ag	Cassette	Secretions	10T
VE90053	Feline Pancreatic Lipase Rapid Test	FPL	Cassette	Serum/Plasma	10T
VE90054	Feline Fecal Occult Blood Rapid Test	F. FOB	Cassette	Feces	10T

Porcine

No.	Product name	Abbr.	Format	Specimen	Specification
VE90055	African Swine Fever Virus Antibody Rapid Test	ASFV Ab	Cassette	Serum/Plasma	20T
VE90056	African Swine Fever Virus Antigen Rapid Test	ASFV Ag	Cassette	Whole Blood/Serum/Plasma	20T
VE90057	Transmissible Gastroenteritis Virus Antigen Rapid Test	TGEV Ag	Cassette	Feces	20T
VE90058	Swine Foot and Mouth Diseases Virus Serotype-O Antibody Rapid Test	S. FMDV-O Ab	Cassette	Serum/Plasma	20T
VE90059	Swine Foot and Mouth Diseases Virus Serotype-A Antibody Rapid Test	S. FMDV-A Ab	Cassette	Serum/Plasma	20T
VE90060	Porcine Epidemic Diarrhea Virus Anti-IgA Rapid Test	PEDV IgA	Cassette	Serum/Plasma/Milk	20T
VE90061	Porcine Circovirus Antibody Rapid Test	PCV Ab	Cassette	Serum/Plasma	20T
VE90062	Porcine Pseudorabies Virus-gE Antibody Rapid Test	PRV-gE Ab	Cassette	Serum/Plasma	20T
VE90063	Porcine Pseudorabies Virus-gB Antibody Rapid Test	PRV-gB Ab	Cassette	Serum/Plasma	20T
VE90064	Porcine Reproductive and Respiratory Syndrom Virus Antibody Rapid Test	PRRSV Ab	Cassette	Serum/Plasma	20T
VE90065	Classical Swine Fever Virus Antibody Rapid Test	CSFV Ab	Cassette	Serum/Plasma	20T

Bovine

No.	Product name	Abbr.	Format	Specimen	Specification
VE90066	Bovine Foot and Mouth Diseases Virus Serotype-0 Antibody Rapid Test	B. FMDV-0 Ab	Cassette	Serum/Plasma	20T
VE90067	Bovine Foot and Mouth Diseases Virus Serotype-A Antibody Rapid Test	B. FMDV-A Ab	Cassette	Serum/Plasma	20T
VE90068	Brucella Antibody Rapid Test	BRU Ab	Cassette	Serum/Plasma	20T
VE90069	Bovine Viral Diarrhea Virus Antibody Rapid Test	BVDV Ab	Cassette	Serum/Plasma	20T
VE90070	Infectious Bovine Rhinotracheitis Antibody Rapid Test	IBR Ab	Cassette	Serum/Plasma	20T

Avian

No.	Product name	Abbr.	Format	Specimen	Specification
VE90071	Newcastle Disease Virus Antigen Rapid Test	NDV Ag	Cassette	Secretions/Cloaca	20T
VE90072	Avian Influenza Virus Antigen Rapid Test	AIV Ag	Cassette	Secretions/Cloaca	20T
VE90073	Avian Influenza Virus H5 Antigen Rapid Test	AIV H5 Ag	Cassette	Secretions/Cloaca	20T
VE90074	Avian Influenza Virus H7 Antigen Rapid Test	AIV H7 Ag	Cassette	Secretions/Cloaca	20T
VE90075	Avian Influenza Virus H9 Antigen Rapid Test	AIV H9 Ag	Cassette	Secretions/Cloaca	20T


Other

No.	Product name	Abbr.	Format	Specimen	Specification
VE90076	Brucella Antibody Rapid Test	BRU Ab	Cassette	Serum/Plasma	20T
VE90077	Clostridium Perfringens Antibody Rapid Test	CLP Ab	Cassette	Serum/Plasma	20T
VE90078	Clostridium Spoilage Antibody Rapid Test	CLS Ab	Cassette	Serum/Plasma	20T
VE90079	Peste des Petits Ruminants Antibody Rapid Test	PPR Ab	Cassette	Serum/Plasma	20T
VE90080	Foot and Mouth Diseases Virus Non-structural Protein 3ABC Antibody Rapid Test	FMDV (3ABC) Ab	Cassette	Serum/Plasma	20T


Test Procedure



- 1** Peel back the seal on the top of the extraction tube containing.
Place the extraction tube in the workstation.



- 2** Use a cotton swab to take a feline fecalor anal swab sample.



- 3** Rotate the swab for at least 10+ times to make sure the specimen is collected effectively and push down the cap onto the extraction tube.



- 4** Add three drops (about 90 ul) of each specimen vertically into the each specimen well (S) of the test device.
Read the result in 10 minutes.

