



# Laboratory Supply Catalog

Volume 24.0



$\text{pH} = \text{pH}[\text{H}^+]$	$[\text{OH}^-]$	Alpha(A)	PH
7.403 68E-08	2.51E-07	0.201	0.201
7.602 51E-08	3.98E-07	0.285	0.285
7.001 00E-08	1.00E-06	0.500	0.500
6.000 98E-09	2.51E-06	0.715	0.715
5.000 158E-09	6.31E-06	0.953	0.953

Plasticware • Glassware • Equipment



## PCR / MICROTUBE RACK & WORKSTATION


Workstation rack holds 5 different sizes of microtubes, screwcap microtubes, PCR tubes and microcentrifuge tubes. Available in blue, green and red.

- Polypropylene (PP)
- Autoclavable at 121°C for 20 minutes

Item #	Color	Unit
111780B	Blue	Each
111780G	Green	Each
111780R	Red	Each

Back row features locking system that holds tubes securely for one-handed opening and closing

Slots and tabs on sides allow multiple racks to be attached to one another



## WORKSTATION RACKS


Designed for use with Globe Scientific sample tubes and CryoClear™ cryogenic vials. A special design allows the self-standing tubes and vials to lock into the rack for easy one-handed opening and closing of caps.

- 40-Place (10 x 4 format)
- Polycarbonate (PC)
- Anti-slip feet
- Autoclavable at 121°C, 15 psig (1 bar) for 20 minutes
- Temperature range: -196°C to +121°C
- Dimensions: 215 x 100 x 25mm (8.5" x 4" x 1")

Item #	Color	Unit
3049B	Blue	Each
3049G	Green	Each
3049R	Red	Each



A broad selection of PCR work racks and cold prep racks are featured on pages 232 & 233.



## FLOAT RACKS FOR 1.5ML AND 2.0 ML MICROTUBES

- Autoclavable polypropylene (PP) construction for strength, buoyancy and temperature resistance
- Racks float with full load of filled tubes
- Molded handle in center of rack for comfort and support during transport
- Legs allow rack to stand upright for use on benchtop
- Use in ice, water baths or on the benchtop
- 65mm round rack fits most 400mL beakers
- Commonly found in cloning, microbiology, and immunology labs

Dimensions (L x W x H):

Square: 95 x 95 x 80mm (3.75" x 3.75" x 3")

Round: 65 dia. x 75mm (2.5 dia. x 3")

Item#	Description	Places	Color	Unit
456395	Square	16 Place	White	4
456398	Round	8 Place	White	4





### DISPENSERS FOR GLASS CULTURE TUBES

These convenient dispensers are made from heavy duty white painted metal. Each holds 1 box of 250 tubes. The front acrylic cover pulls up to dispense the tubes.

- Conserve valuable benchtop real estate
- Keep tubes clean until point of use
- Safely remove only what you need

Item #	Description	Unit
1590	For 10 x 75mm and 12 x 75mm tubes	Each
1592	For 13 x 100mm tubes	Each

1590, 1592 — Work great with Globe Scientific glass tubes! See page 278.



### SEROLOGICAL PIPETTE STORAGE RACK

- Four angled compartments protect and organize glass or plastic serological pipets
- Durable ABS construction for scratch resistance and easy cleaning
- Rack can stand on flat surface without mounting
- Powerful magnets on both side for secure mounting to most metal surfaces
- Two included accessory attachments keep pipet fillers and bulbs within reach
- Side windows allow for easing viewing of pipet size and inventory

Dimensions (L x W x H): 290 x 85 x 395mm (11.4" x 3.3" x 15.5")

Item#	Description	Color	Unit
1685B	Pipette storage rack ABS	Blue	Each
1685R	Pipette storage rack ABS	Red	Each
1685W	Pipette storage rack ABS	White	Each



### HEAVY-DUTY WALL MOUNT DRYING/DRAINING RACK

This 72 place laboratory drying rack is perfect for draining and drying of all types of labware. The rack comes with a mounting kit to hang on the wall. The back plate of the support is made from a single mold of 4mm rigid high impact polystyrene (HIPS). The peg sockets have closed ends eliminating leakage and consequent biohazards. The seamless drip channel collects waste which drains via the supplied drain tube.

The 72 (push in) pegs can accommodate items with a neck bore in excess of 15mm. They are easily removable for cleaning or to accommodate odd shaped items. Also included are 11 small pegs to accommodate small diameter items, such as test tubes. This durable rack will not rust or stain.

- 72 interchangeable pegs
- Mounting kit and template for hanging
- Drip channel with draining hose
- Hooks for joining together additional racks

Dimensions: 24.75" high x 17.75" wide x 4.5" deep

Pegs: 3.7" (L) x 0.6" (D) (95mm x 15mm)

Small Pegs: 3.7" (L) x 0.024" (D) (95mm x 6mm)

Item #	Description	Unit
600213	Drying rack, 72-place	Each
601213	Additional small pegs for drying rack, 11/pack	11

