



Inorganic Disinfection #2, WS

One 24 mL screw-cap vial yields up to 4 liters after dilution.



VARIANTS

Item no.	Name	Product Type	Study Frequency	Study Type	Accreditations
5262	Inorganic Disinfection 2 QC	CRM		Water Supply (WS)	ISO 17034

Item no.	Name	Product Type	Study Frequency	Study Type	Accreditations
Specifications					
Product Family		Environmental			
Product Category		Water Supply			
Product Sub Category		Inorganics			
Methods		ASTM D6581-00(05), ASTM D6919-03, EPA 300.0, EPA 300.1, EPA 317.0, EPA 320.1, EPA 326, EPA 9211, SM 4110 B, SM 4500-Br ⁻ D			
Storage		20 - 25C (Refrigerator)			
Shelf Life from date of manufacture (months)		30			
Minimum shelf life on ship date (months)		12			
Container Type		24mL Amber Glass Vial			
Components/Pkg		1			
Volume		24 mL			
Weight		0.05 lbs			
Application		Environmental			
5260	Inorganic Disinfection 2	PT	Monthly	Water Supply (WS)	ISO 17043

Item no.	Name	Product Type	Study Frequency	Study Type	Accreditations
Specifications					
Product Family		Environmental			
Product Category		Water Supply			
Product Sub Category		Inorganics			
Methods		ASTM D6581-00(05), ASTM D6919-03, EPA 300.0, EPA 300.1, EPA 317.0, EPA 320.1, EPA 326, EPA 9211, SM 4110 B, SM 4500-Br ⁻ D			
Storage		20 - 25C (Refrigerator)			
Shelf Life from date of manufacture (months)		30			
Minimum shelf life on ship date (months)		12			
Container Type		24mL Amber Glass Vial			
Components/Pkg		1			
Volume		24 mL			
Weight		0.05 lbs			
Application		Environmental			
5262QR	Inorganic Disinfection #2, WS QR	QR		Water Supply (WS)	ISO 17043

Item no.	Name	Product Type	Study Frequency	Study Type	Accreditations
Specifications					
	Product Family	Environmental			
	Product Category	Water Supply			
	Product Sub Category	Inorganics			
	Methods	ASTM D6581-00(05), ASTM D6919-03, EPA 300.0, EPA 300.1, EPA 317.0, EPA 320.1, EPA 326, EPA 9211, SM 4110 B, SM 4500-Br ⁻ D			
	Storage	20 - 25C (Refrigerator)			
	Shelf Life from date of manufacture (months)	30			
	Minimum shelf life on ship date (months)	12			
	Container Type	24mL Amber Glass Vial			
	Components/Pkg	1			
	Volume	24 mL			
	Weight	0.05 lbs			
	Application	Environmental			

ANALYTES

Analyte	Low Conc	High Conc	Units
Bromate	7	50	µg/L
Bromide	50	300	µg/L