



Solids Concentrate, WS

One 24 mL screw-cap vial with a powder yields 1 liter after dilution.




VARIANTS

| Item no. | Name | Product Type | Study Frequency | Study Type | Accreditations |
|----------|-----------------------------------|--------------|-----------------|-------------------|----------------|
| 5152 | Solids, Concentrate, PotableWatR™ | CRM | | Water Supply (WS) | ISO 17034 |

Questions? We are happy to help.




| Item no. | Name | Product Type | Study Frequency | Study Type | Accreditations |
|--|-------------------------|--|-----------------|-------------------|----------------|
| Specifications | | | | | |
| Product Family | | Environmental | | | |
| Product Category | | Water Supply | | | |
| Product Sub Category | | Inorganics | | | |
| Methods | | EPA 160.1, EPA 160.2, EPA 160.3, SM 2540 B, SM 2540 C, SM 2540 D, USGS I-3765-85 | | | |
| Storage | | 20 - 25C (Refrigerator) | | | |
| Shelf Life from date of manufacture (months) | | 42 | | | |
| Minimum shelf life on ship date (months) | | 12 | | | |
| Container Type | | 24mL Clear Glass Vial | | | |
| Components/Pkg | | 1 | | | |
| Volume | | 24 mL | | | |
| Weight | | 0.05 lbs | | | |
| Application | | Environmental | | | |
| 5150 | Solids, Concentrate, WS | PT | Monthly | Water Supply (WS) | ISO 17043 |



Questions? We are happy to help.



| Item no. | Name | Product Type | Study Frequency | Study Type | Accreditations |
|--|---------------------------|--|-----------------|-------------------|----------------|
| Specifications | | | | | |
| Product Family | | Environmental | | | |
| Product Category | | Water Supply | | | |
| Product Sub Category | | Inorganics | | | |
| Methods | | EPA 160.1, EPA 160.2, EPA 160.3, SM 2540 B, SM 2540 C, SM 2540 D, USGS I-3765-85 | | | |
| Storage | | 20 - 25C (Refrigerator) | | | |
| Shelf Life from date of manufacture (months) | | 42 | | | |
| Minimum shelf life on ship date (months) | | 12 | | | |
| Container Type | | 24mL Clear Glass Vial | | | |
| Components/Pkg | | 1 | | | |
| Volume | | 24 mL | | | |
| Weight | | 0.05 lbs | | | |
| Application | | Environmental | | | |
| 5152QR | Solids Concentrate, WS QR | QR | | Water Supply (WS) | ISO 17043 |



Questions? We are happy to help.



| Item no. | Name | Product Type | Study Frequency | Study Type | Accreditations |
|----------------|--|--|-----------------|------------|----------------|
| Specifications | | | | | |
| | Product Family | Environmental | | | |
| | Product Category | Water Supply | | | |
| | Product Sub Category | Inorganics | | | |
| | Methods | EPA 160.1, EPA 160.2, EPA 160.3, SM 2540 B, SM 2540 C, SM 2540 D, USGS I-3765-85 | | | |
| | Storage | 20 - 25C (Refrigerator) | | | |
| | Shelf Life from date of manufacture (months) | 42 | | | |
| | Minimum shelf life on ship date (months) | 12 | | | |
| | Container Type | 24mL Clear Glass Vial | | | |
| | Components/Pkg | 1 | | | |
| | Volume | 24 mL | | | |
| | Weight | 0.05 lbs | | | |
| | Application | Environmental | | | |

ANALYTES

| Analyte | Low Conc | High Conc | Units |
|---------------------------------|----------|-----------|-------|
| Non-Filterable Residue (TSS) | 23 | 100 | mg/L |
| Total Dissolved Solids at 180°C | 100 | 1000 | mg/L |
| Total Solids at 105°C | 123 | 1100 | mg/L |

×

Questions? We are happy to help.

