



## **Product Specification Sheet**

### **TS/19-G250**

Polystyrene shaft swab with polyester tip sealed in tube.

Product Code	TS/19-G250
Description	Polystyrene shaft swab with polyester tip sealed in tube.
Intended Use	Culture swab; Patient sampling & general Laboratory use
Components	Polypropylene Tube, Blue Standard Cap, Polystyrene/Polyester tip swab, Non-toxic glue, Label.
Contents	Not applicable.
Method of Sterilisation	Sterilised by gamma irradiation at dosage of 25-53.4kGy
Shelf Life	24 months from date of manufacture
Storage Conditions	4-25°C. Store in dry and ventilated area; protect from light, heat and moisture in tightly closed original packaging
Precautions and Limitations	Not applicable.
Pack Size	250 in a box
Other	Quality control testing details on final product: i) Microbial efficacy test using indicator organism(s) ii) Validation of irradiation sterility.

This specification (and formulation where applicable) are intended for guidance only and may be adjusted to meet product defined criteria.

You should always refer to a professional with experience in your intended application if in any doubt.

Issue #: 3 Date: 24.04.2020  
ProductSpec\_TS19-G250