

NEW!

TransferTop™

Urine Collection & Direct Transfer

3856L — 4oz (120mL)
TransferTop™ Collection and Transfer Cup with integrated transfer mechanism (graduated to 100mL)

TransferTop™ Urine Collection and Direct Transfer Cups feature an integrated transfer mechanism for “closed system” transfer of urine samples directly from the cup to a vacuum sample tube. With the TransferTop™ System, there is no need for laboratory personnel to open the collection cup and pour the sample into a tube for processing. Instead, the safety label is peeled back, and the vacuum tube inserted for simple, safe, efficient, and reliable sample transfer.

TransferTop™ Closed System for Urine Transfer:

- Requires no open transferring of samples for reduced spills, improved lab hygiene, and limited exposure to dangerous samples
- Reduces risk of contamination for improved test results
- Directly transfers samples to a tube suitable for analysis and transport
- Works with most vacuum urine tubes
- Provides consistent sample transfer volumes when used with like vacuum tubes
- Improves efficiencies and sample throughput when compared to traditional open systems

Item #	Description	Unit
3852	Urine Collection & Direct Transfer Cup, 2oz (60mL) capacity, STERILE, individually wrapped	500
3856L	Urine Collection & Direct Transfer Cup, 4oz (120mL) capacity, STERILE, bulk	300
3841	Urine Transfer Straw, 3" (7.5cm) straw, non-sterile, 100/bag, 8 bags/case	800
108010	3 Liter, 24 hour urine collection container with attached TransferTop™ screwcap, amber	30



3852 — 2oz (60mL) TransferTop™ Collection and Transfer Cup



3841 — Urine Transfer Straw for closed transfer of urine when drawing from a standard collection container



108010 — 3L TransferTop™ 24 Hour Urine Collection Container features an integrated mechanism for “closed system” transfer of urine directly from the container to vacuum tubes.



ALL NEW!



24 HOUR URINE COLLECTION CONTAINERS

New 24-hour urine collection containers offer sizes, configurations, and features required to meet most 24-hour collection needs. Upright configurations for efficient storage and low form designs for stability are all designed with easy to read, molded graduations and leak resistant closures. Amber polyethylene construction protects samples from light and molded handles make handling easy. The new TransferTop™ option includes an integrated mechanism for “closed system” transfer to vacuum urine sample tubes.

- Amber polyethylene (PE) construction
- Easy to read molded graduations
- Molded handles for easy handling
- Leak resistant closures
- TransferTop™ option for improved lab safety and hygiene

Item #	Description	Unit
108010	3 Liter, with attached TransferTop™ screwcap, amber	30
108020A	2 Liter, with attached screwcap, amber	54
108025A	2.5 Liter, with attached screwcap, amber	30
108030A	3 Liter, with attached screwcap, amber	30

108010— 3L TransferTop™ 24 Hour Urine Collection Container features an integrated mechanism for “closed system” transfer of urine directly from the container to vacuum tubes.

diamond ESSENTIALS

Quality. Value. No Compromises.

NEW!



SPECIMEN CONTAINERS FOR ON-SITE COLLECTION AND TESTING

4 oz (120mL) Urine Collection Container with Full Turn Screw Cap

These 4 ounce urine collection containers are ideal for use in physician's offices and hospital settings. The containers feature a full turn screw cap for better leak resistance. They have molded graduation marks and a frosted writing area. The individually wrapped sterile containers also come complete with an affixed patient ID label.

- Graduated to 120mL
- Frosted writing area
- Clarified polypropylene (PP)

Sterile & Aseptic Containers:

- Individually wrapped
- Trilingual patient ID label (English, Spanish & French)

Material: Container = Polypropylene (PP) • Cap = Polyethylene (PE)

Dimensions: Container: 57.6 x 72.0mm • Cap: 62mm O.D.

Item #	Description	Unit
5901G	Full turn green attached screw cap, ID label, individually wrapped, STERILE	250
5901R	Full turn red attached screw cap, ID label, individually wrapped, STERILE	250
5902G	Full turn green attached screw cap, ID label, individually wrapped, CLEANROOM ISO8 ASEPTIC	250
5902R	Full turn red attached screw cap, ID label, individually wrapped, CLEANROOM ISO8 ASEPTIC	250
5905G	Full turn green separate screw cap, non-sterile	500
5905R	Full turn red separate screw cap, non-sterile	500



5912



5913



Tite-Rite™ leak resistant containers with patient I.D. and sterility assurance labels



SPECIMEN CONTAINERS AND COLLECTION CUPS

4 oz (120mL) Urine Collection Container with 1/4 Turn Screw Cap

Ideal for use in physician's offices and hospital settings. The containers have molded graduation marks and a frosted writing area.

Sterile Containers:

- Graduated to 100mL
- Individually wrapped
- Trilingual patient ID label (English, Spanish & French)

Non-Sterile Containers:

- Graduated to 100mL
- Frosted writing area

Material: Container = Polypropylene (PP) • Cap = Polyethylene (PE)

Dimensions: Container: 56.8 x 72.5mm • Cap: 61mm O.D.

Item #	Description	Unit
5912	Red screw cap, ID label, individually wrapped, ASEPTIC	100
5913	Green screw cap, ID label, individually wrapped, ASEPTIC	100
5914	Red screw cap, non-sterile	100
5915	Red screw cap, non-sterile	500
5915G	Green screw cap, non-sterile	500
5916	Red screw cap, affixed thermometer strip, non-sterile	500

LEAK RESISTANT CONTAINERS

- Ideal for hospital pneumatic tube systems
- Tite-Rite caps are produced from polyethylene and contain an inner liner designed to compress against the uppermost walls of the container, providing an air-tight seal

Container: Polypropylene (PP) • Cap: Polyethylene (PE)

For Urine Collection and Transport

- Unbroken seal guarantees sterility of interior at point of use


Item #	Description	Unit
6220	90mL (48mm) container, with sterility assurance label, STERILE	400
6526	90mL (53mm) container, with sterility assurance label, STERILE	300
6528	120mL (53mm) container, with sterility assurance label, STERILE	300
6529	120mL (53mm) container, with sterility assurance label, individually wrapped, STERILE	100
6219	90mL (48mm) container with label, non-sterile	400

LEAK-RESISTANT DRUG TESTING CONTAINERS

60mL and 90mL Container with Thermometer Strip

These urine specimen containers come with attached thermometer strip, tamper evident ID label and white screw cap. The 90mL containers are available in 48mm regular and 53mm wide mouth configurations.

Item #	Description	Unit
6523TP	60mL (45mm) container, STERILE, 400/case	500
6238	90mL (48mm) container, STERILE, 400/case	400
6236	90mL (53mm) wide mouth container, STERILE, 300/case	300



3800 — Quick-Read Precision Grid urinalysis slide

MICROSCOPIC URINALYSIS

Quick-Read™ Precision Grid Urinalysis Microscope Slide


The Quick-Read™ Precision Grid slide is an innovative product designed to provide accuracy, uniformity and safety in the microscopic examination of urinary sediment. Made of optically clear acrylic for optimal viewing, the slide consists of ten chambers. Each chamber has a grid with 81 well defined, pre-measured squares. Each pre-measured square holds a specific volume of urinary sediment which facilitates the convenient and rapid microscopic examination of the cellular elements in each specimen.

Item #	Description	Unit
3800	Quick-Read™ Precision Grid microscope slide, 100 slides/box for 1000 tests	100

Quick-Read™ Precision Grid Urinalysis System Kits

These convenient kits include all necessary components for the microscopic examination of urinary sediment. Each kit contains supplies for 1000 tests.

Item #	Description	Unit
3817	Quick-Read Precision Grid slide, Quick-Prep tube & pipette kit, plug cap & stain (for 1000 tests)	Case
3819	Quick-Read Precision Grid slide, Quick-Prep tube & pipette kit, plug cap (for 1000 tests)	Case
3818	Quick-Read Precision Grid slide, Quick-Prep tube & pipette kit (for 1000 tests)	Case





3805 — Quick-Read Precision Cell urinalysis microscope slide

MICROSCOPIC URINALYSIS

Quick-Read™ Precision Cell Urinalysis Microscope Slide

This innovative slide has been designed to provide accuracy, uniformity and safety in the microscopic examination of urinary sediment.



Made from optically clear acrylic, the slide features 10 chambers with integrated coverslips. Within each chamber there are 18 circles that are all together visible at 100x magnification. These premeasured circles contain specific volumes of urinary sediment within its circumference that allows for the easy calculation of cellular elements present in 1mL of urine.

Item #	Description	Unit
3805	Quick-Read™ microscope slide, 100 slides/box for 1000 tests	100

Quick-Read™ Precision Cell Urinalysis System Kits

All necessary components for the microscopic examination of urinary sediment. Each kit contains supplies for 1000 tests.

Item #	Description	Unit
3827	Quick-Read Precision Cell slide, Quick-Prep tube & pipette kit, plug cap & stain (for 1000 tests)	Case
3829	Quick-Read Precision Cell slide, Quick-Prep tube & pipette kit, plug cap (for 1000 tests)	Case
3828	Quick-Read Precision Cell slide, Quick-Prep tube & pipette kit (for 1000 tests)	Case



Area: 63mm²
Volume: 6.3μL



Total Area: 1.0mm²
Total Volume: 0.1μL



Area: 0.111mm²
Volume: 0.011μL

Quick-Prep™ Tube & Pipette Kit

This tube and pipette combination has been specially designed for use with the Quick-Read microscope slides. The 12mL polystyrene (PS) centrifuge tube is graduated at the 10mL mark for convenient filling of urine without prior measurement. The pipette has a specially designed bulb used to hold 1mL of urine sediment in place while the remaining urine is decanted from the tube.


Item #	Description	Unit
3825	Quick-Prep™ tube & pipette kit, 500/pack, 2 packs/case (for 1000 tests)	1000

Accessories

Item #	Description	Unit
3810	Sternheimer Malbin Urine Sediment Stain, 50mL/bottle	Each
456106	10-Place rack, holds 10 Quick-Prep tubes	Each
112010	Quick-Prep 12mL (16 x 100mm) centrifuge tube, PS, 500/pack, 4 packs/case	2000
116142	16mm plug cap for #112010 centrifuge tube, white, 1000/pack	1000
3826	Quick-Prep transfer pipette with sediment bulb, 1mL draw, 500/pack, 6 packs/case	3000



3825 — Quick-Prep tube and pipette



URI-PAK™ URINE COLLECTION KITS

The Uri-Pak™ Urine Collection Kit offers all the components necessary for specimen collection, handling, transport and testing. Uri-Pak is used in hospitals, laboratories, physicians' offices and nursing homes throughout the world.

The Uri-Pak System is economically priced and may be purchased as a complete kit or by its separate components. Custom kits and private label packaging are available.

The kit consists of five plastic bags each containing:

- 1 100 — flared-top polystyrene centrifuge tubes
- 2 100 — polyethylene snap-caps
- 3 100 — adhesive patient ID labels
- 4 100 — 3.5 oz. collection cups
- 5 Optional 10-place disposable test tube rack

Custom kits are available. Call to learn how we can customize a solution for your particular needs.



112018R— Uri-Pak® Kit with 12mL Ames®-type tubes

Item #	Description	Unit
112018R	Uri-Pak Kit with 12mL Ames®-type tube (#112017), blue caps, cups, labels & rack	500
112018	Same as above but without rack	500
112019R	Uri-Pak Kit with 12mL Ames®-type tube (#112017), yellow caps, cups, labels & rack	500
112019	Same as above but without rack	500



112015
Round bottom



112017
Ames®-type

Tubes shown actual size



112030
Whale®-type



112010/112009

URINE CENTRIFUGE TUBES

Flared-Top 12 mL Centrifuge Tubes

This line of centrifuge tubes has been designed specifically for urinalysis procedures. The round bottom, Ames® and Whale® type tubes feature a flared opening which is designed for easy pour-off of liquids.

- Polystyrene (PS)
- Spin rate: 1500 RCF
- Dimensions: 20 x 115mm

Item #	Description	Unit
112015	12mL urine tube, round bottom, graduated, flared opening, 2 x 1000/cs	2000
112017	12mL urine tube, (Ames -Type), graduated, flared opening, 3 x 500/cs	1500
112030	12mL urine tube (Whale -Type), w/sediment bulb, graduated, flared opening, 3 x 500/cs	1500
113136	Snap-cap, blue, 3 x 500/case	1500
113137	Snap-cap with sanitary grip, yellow, 3 x 500/case	1500

12 mL Centrifuge Tube

- Dimensions: 16 x 100mm
- Graduated at 0.5, 1, 2.5, 5 and 10mL
- Conical bottom
- Polystyrene (PS)
- Spin rate: 1500 RCF

Item #	Description	Unit(s)
112010	Polystyrene (PS), 500/bag, 4 bags/case	500 & 2000
112009	Polystyrene (PS), with separate plug cap	1000
116142	16mm plug stopper, white	1000
116152C	16mm flanged plug stopper, natural (additional colors available)	1000

Caps for 12mL flared-top centrifuge tubes:



113136



113137

Caps for 12mL centrifuge tube:



116142



116152C

12 mL Centrifuge Tubes — Graduated to 12 mL

- Dimensions: 16.5 x 107.5mm
- Graduated at 0.25, 0.5, 1, 2, 2.5, 3, 4, 5, 6, 7, 8, 9, 10, 11 and 12mL
- Conical bottom
- Polystyrene (PS)
- Spin rate: 4000 x g

Item #	Description	Unit
112050	Polystyrene (PS), graduated, 100/bag	1000
6263C	Plug cap, natural	1000
6263B	Plug cap, blue	1000
6263G	Plug cap, green	1000
6263R	Plug cap, red	1000
6263Y	Plug cap, yellow	1000
6263W	Plug cap, white	1000

112050
Polystyrene (PS)



Tube shown actual size



6263C

15 mL Centrifuge Tube

- Dimensions: 16.6 x 120mm
- Graduated from 0.5mL to 15mL in 0.5mL increments
- Conical bottom
- Polystyrene (PS)
- Spin rate: 1500 RCF

Item #	Description	Unit(s)
6260	Polystyrene (PS), graduated, 250/pack, 4 packs/case	250 & 1000
6265	Polystyrene (PS), graduated, 100/dispenser box, 10 boxes/case	1000
6262	Polystyrene (PS), graduated, with separate plug cap included	1000
6263C	Plug cap for 15mL tubes, natural	1000

6265 — Available in convenient benchtop dispenser boxes.



6260/6265/6262



Tube shown actual size



6263C



Contáctenos:
310 6888259 • 314 2852622
www.labcarecolombia.com
ventas@labcarecolombia.com