

diamond ESSENTIALS

Quality. Value. No Compromises.

NEW!



SPECIMEN CONTAINERS FOR ON-SITE COLLECTION AND TESTING

4 oz (120mL) Urine Collection Container with Full Turn Screw Cap

These 4 ounce urine collection containers are ideal for use in physician's offices and hospital settings. The containers feature a full turn screw cap for better leak resistance. They have molded graduation marks and a frosted writing area. The individually wrapped sterile containers also come complete with an affixed patient ID label.

- Graduated to 120mL
- Frosted writing area
- Clarified polypropylene (PP)

Sterile & Aseptic Containers:

- Individually wrapped
- Trilingual patient ID label (English, Spanish & French)

Material: Container = Polypropylene (PP) • Cap = Polyethylene (PE)

Dimensions: Container: 57.6 x 72.0mm • Cap: 62mm O.D.

Item #	Description	Unit
5901G	Full turn green attached screw cap, ID label, individually wrapped, STERILE	250
5901R	Full turn red attached screw cap, ID label, individually wrapped, STERILE	250
5902G	Full turn green attached screw cap, ID label, individually wrapped, CLEANROOM ISO8 ASEPTIC	250
5902R	Full turn red attached screw cap, ID label, individually wrapped, CLEANROOM ISO8 ASEPTIC	250
5905G	Full turn green separate screw cap, non-sterile	500
5905R	Full turn red separate screw cap, non-sterile	500

SPECIMEN CONTAINERS FOR ON-SITE COLLECTION AND TESTING

4 oz (120mL) Urine Collection Container with 1/4 Turn Screw Cap

Our popular 4 ounce urine collection container is ideal for use in physician's offices and hospital settings. The containers have molded graduation marks and a frosted writing area.

Sterile Containers:

- Graduated to 100mL
- Individually wrapped
- Trilingual patient ID label (English, Spanish & French)

Non-Sterile Containers:

- Graduated to 100mL
- Frosted writing area

Material: Container = Polypropylene (PP) • Cap = Polyethylene (PE)

Dimensions: Container: 56.8 x 72.5mm • Cap: 61mm O.D.

Item #	Description	Unit
5912	Red screw cap, ID label, individually wrapped, ASEPTIC	100
5913	Green screw cap, ID label, individually wrapped, ASEPTIC	100
5914	Red screw cap, non-sterile	100
5915	Red screw cap, non-sterile	500
5915G	Green screw cap, non-sterile	500
5916	Red screw cap, affixed thermometer strip, non-sterile	500



5912




5913

5 oz Wide Mouth Container with Snap Cap

- Extra wide 70mm (2.75") top inner diameter
- Height with cap: 67.7mm, without cap: 66.4mm
- Molded graduations in both ounces and milliliters

Material: Container = Polypropylene (PP) • Cap = Polyethylene (PE)

Item #	Description	Unit
5917	With snap cap, 300/case	300



5917

3.5 oz Collection Cup (No Lid)

Polystyrene (PS)

Item #	Description	Unit
6207	Plastic flexible cup, 100/pk, 25 packs/case	2500



6207



Labcare
de Colombia
www.labcarecolombia.com





LEAK RESISTANT DRUG TESTING CONTAINERS

60mL and 90mL Containers with Thermometer Strip

These urine specimen containers come with attached thermometer strip, tamper evident ID label and white screw cap. The 90mL containers are available in 48mm regular and 53mm wide mouth configurations.

Item #	Description	Unit
6523TP	60mL (45mm) container, STERILE, 500/case	500
6238	90mL (48mm) container, STERILE, 400/case	400
6236	90mL (53mm) wide mouth container, STERILE, 300/case	300

Both sizes shown with attached tamper evident ID labels and thermometer strips



6.5 oz. collection cups (5920-5929)

6.5 OZ (190ML) COLLECTION CUP

- Graduated
- Pour Spout
- Packaged 25 per sleeve, 20 sleeves per case

Item #	Description	Unit
5920	Polypropylene (PP) beaker, 25/pk, 20 pks/case	500
5924	Polypropylene (PP) beaker, with paper lids, 25/pk, 20 pks/case	500
5926	Polypropylene (PP) beaker, with attached thermometer strip	500
5921	Polystyrene (PS) beaker, 25/pk, 20 pks/case	500
5925	Polystyrene (PS) beaker, with paper lids, 25/pk, 20 pks/case	500
5929	Paper lid for 6.5 oz beakers, 25/pk, 20 pks/case	500



112215

THERMOMETER STRIPS - LIQUID CRYSTAL

The self-adhesive liquid crystal thermometer strips are ideal for use with all drug testing collection containers. They are graduated in both Fahrenheit (90° - 100°F) and Celsius (32°-38°C).

Item #	Description	Unit
112215	Liquid crystal thermometer strips, 4000 strips per roll	4000

PATIENT IDENTIFICATION LABELS

Self-adhesive patient identification labels for use on all specimen collection containers. The labels come 25 per sheet, 80 sheets per case.

Item #	Description	Unit
112192	Patient identification label, self-adhesive, 25/sheet	2000



112192



LEAK RESISTANT CONTAINERS

- Tite-Rite™ containers are the ultimate in leak-tight storage and transport of liquids
- Ideal for hospital pneumatic tube systems
- Tite-Rite caps are produced from polyethylene and contain an inner liner designed to compress against the uppermost walls of the container, providing an air-tight seal
- The containers are produced from polypropylene and have molded graduations

Container: Polypropylene (PP) • Cap: Polyethylene (PE)

Non-Sterile, Bulk Packaging

Item #	Size	Diameter	Height	Unit
6518	20mL (0.66 oz)	35mm	44mm	1000
6520	40mL (1.33 oz)	48mm	41mm	600
6522	60mL (2 oz)	48mm	60mm	500
6524	90mL (3 oz)	48mm	76mm	400
6525	90mL (3 oz)	53mm	63mm	300
6527	120mL (4 oz)	53mm	78mm	300



Tite-Rite™ leak resistant containers

For Urine Collection and Transport — **STERILE with Sterility Assurance Label**

- Unbroken seal guarantees sterility of interior at point of use

Item #	Description	Unit
6220	90mL (48mm) container, with sterility assurance label, STERILE	400
6523	60mL container with label, STERILE	500
6526	90mL (53mm) container, with sterility assurance label, STERILE	300
6528	120mL (53mm) container, with sterility assurance label, STERILE	300
6529	120mL (53mm) container, with sterility assurance label, individually wrapped, STERILE	100



For Urine Collection and Transport — **Non-Sterile**

Item #	Description	Unit
6219	90mL (48mm) container with label, non-sterile	400

Don't forget to see our complete line of specimen transport bags. Several styles and multiple sizes are available.
Pages 24-25.



120mL container with patient I.D. and sterility assurance label, individually wrapped



**Labcare
de Colombia**
www.labcarecolombia.com





109217 — polystyrene (PS), 109220 — polypropylene (PP)

30mL TRANSPORT TUBES

These 30mL screw cap specimen transport tubes are available in transparent polystyrene (PS) or translucent polypropylene (PP). The self-standing tubes feature a conical bottom with skirt.

30mL Self-Standing with Screw Cap

- Dimensions of tube: 25 x 90mm • Dimensions with cap: 30 x 94mm
- Conical bottom with skirt
- Non-graduated
- Non-sterile

Item #	Description	Unit
109217	30mL polystyrene (PS) tube with polypropylene (PP) cap	500
109220	30mL polypropylene (PP) tube with polyethylene (PE) cap	500



109117 — polystyrene (PS), 109120 — polypropylene (PP)

FECAL COLLECTION CONTAINERS

30mL Self-Standing with Attached Spoon

Globe's 30mL screw cap container with spoon is ideal for the collection and transport of fecal specimens. Available in transparent polystyrene or opaque polypropylene, these stool containers come completely assembled and are ready to use.

Cap: Polypropylene • Spoon: Polystyrene

Item #	Description	Unit
109117	30mL, polystyrene, with spoon and screwcap	500
109120	30mL, polypropylene, with spoon and screwcap	500



109224

20mL with Attached Spoon

20mL white stool collection container with blue lid and attached spoon.

Cap and spoon: Polypropylene

Item #	Description	Unit
109224	20mL, polypropylene, 100/bag, 12 bags/case	1200



Labcare
de Colombia
www.labcarecolombia.com



MULTI-PURPOSE CONTAINERS WITH SNAP ON LID**Translucent (To View Contents) with Lid**

Container: High density polyethylene (HDPE)

Snap Cap: Low density polyethylene (LDPE)

Item #	Size	Diameter x Height	Unit
271004*	4 ounce (120mL)	2.95" x 2.23"	300
271008	8 ounce short (240mL)	3.80" x 2.23"	100
271009	8 ounce tall (240mL)	3.25" x 3.48"	100
271016	16 ounce (480mL)	3.90" x 4.36"	100
271032	32 ounce (1000mL)	4.60" x 5.40"	100
271064	64 ounce (2000mL)	5.90" x 6.25"	50
271086	86 ounce (2500mL)	5.90" x 7.60"	25
271172	172 ounce (5000mL)	8.60" x 7.00"	10

*NOTE: The 4 oz container (item #271004) is made from Polypropylene (PP)

Opaque White with Lid

Container: High density polyethylene (HDPE)

Snap Cap: Low density polyethylene (LDPE)

Item #	Size	Diameter x Height	Unit
271008W	8 ounce short (240mL)	3.80" x 2.23"	100
271009W	8 ounce tall (240mL)	3.25" x 3.48"	100
271016W	16 ounce (480mL)	3.90" x 4.36"	100
271032W	32 ounce (1000mL)	4.60" x 5.40"	100
271064W	64 ounce (2000mL)	5.90" x 6.25"	50
271086W	86 ounce (2500mL)	5.90" x 7.60"	25
271172W	172 ounce (5000mL)	8.60" x 7.00"	10

Opaque White with Lid - Extra Heavy-Duty

Item #	Description	Dia. x Height	Unit
270008	8 ounce (240mL), polypropylene (PP)	3.75" x 2.25"	250
270016	16 ounce (480mL), polypropylene (PP)	4" x 3"	250
270032	32 ounce (1000mL), polyethylene (PE)	4.5" x 4.75"	250
270064	64 ounce (2000mL), polyethylene (PE)	6.5" x 5"	150
270086	86 ounce (2500mL), polypropylene (PP)	6" x 7"	90
270172	172 ounce (5000mL), polyethylene (PE)	8.25" x 7"	40
270190	190 ounce (5700mL), polyethylene (PE)	8.25" x 7.5"	20

**HISTOLOGY CONTAINERS —
EXTRA LARGE CAPACITY WITH SCREW CAP**

These extra large capacity containers are ideal for histology specimens. They are bulk packaged with separate yellow screwcaps.

- Non-sterile
- Stackable
- Ribbed cap for easy opening and closing

Container: Polypropylene (PP) • Cap: Polyethylene (PE)

Item #	Size	Diameter x Height	Unit
6540	250mL	3.5" x 2.25"	100
6542	500mL	3.5" x 4.375"	100
6545	1000mL	4.375" x 5.375"	100



Translucent containers — 271004 – 271172



Opaque containers — 271008W – 271172W



**Labcare
de Colombia**
www.labcarecolombia.com

